

Der Procyon-Kompressorvernebler bietet eine schnelle, effiziente und zuverlässige Vernebelungsinhalationstherapie, die für den Heimgebrauch konzipiert ist. Dies ist ein leiser und leistungsstarker Kompressorvernebler, der sich zum Vernebeln der meisten Aerosolmedikamente eignet, die gegen häufige Atemwegserkrankungen verschrieben werden. Er bietet höhere Literdurchflüsse, kürzere Behandlungszeiten und niedrigere Geräuschpegel als viele andere Economy-Kompressorvernebler.

**Leistungs- und Designmerkmale:**

- Zuverlässiger Druck und Durchfluss
- Ein wiederverwendbarer Vernebler kann enthalten sein
- Kleiner Fußabdruck
- Klares, kompaktes Design
- Praktischer Tragegriff
- Mechanische Teile aus Metall

Weight: 2.18kg



Serie/Modell-Nr.	Procyon-Serie / 70802
Gerätemessungen	191 mm * 173 mm * 102 mm
Gewicht	2.18 kg
Elektrische Anforderungen	100-240 VAC, 50-60 Hz
Betriebseräuschpegel	≤ 60 dBA
Maximaler Druck	≤ 35 psi
Freier Luftstrom des Kompressors	≤ 6 l/min
Partikelgröße	2µm - 5µm
Verneblertrate	≥ 0.2 ml/min
Restflüssigkeit	<0.5 ml




## PORTABLE COMPRESSOR NEBULIZER

Penetrate into Bronchus & Lung 2-5µm



## MICRON NEBULIZATION FINE PARTICLES



<p><b>3.6-4.0µm</b> Median particle size MMAD</p>	<p><b>75%</b> Percentage of effective particles</p>
---	---

## QUIET ENVIRONMENT SILENT ATOMIZATION

Low noise design, easy for children to sleep



Nebulizer is suitable for auxiliary treatment of most **respiratory diseases**



1. Nebulzier machine
2. Medicine cup
3. Air tube
4. Adult mask
5. Child mask
6. Mouthpiece
7. Elastic strap

