

El nebulizador compresor Procyon proporciona una terapia de inhalación nebulizada rápida, eficiente y confiable, diseñada para uso doméstico. Se trata de un nebulizador compresor silencioso y resistente, adecuado para nebulizar la mayoría de los medicamentos en aerosol recetados para enfermedades respiratorias comunes. Ofrece mayores flujos de litros, tiempos de tratamiento más cortos y niveles de sonido más bajos que muchos otros nebulizadores de compresor económicos.

Características de rendimiento y diseño:

- Presión y flujo confiables
- Se puede incluir un nebulizador reutilizable
- Pequeña huella de pie
- Diseño pequeño y compacto
- Cómoda asa de transporte
- Piezas mecánicas metálicas.

Weight: 2.18kg



Serie / Modelo No.	Serie Procyon / 70802
Dimensiones de la unidad	191mm*173mm*102mm
Peso	2.18 kg
Requisitos eléctricos	100-240 VCA, 50-60 HZ
Nivel de ruido de funcionamiento	≤ 60 dBA
Presión máxima	≤ 35 psi
Flujo de aire libre del compresor	≥ 6l/min
Tamaño de partícula	2 μm - 5 μm
Tasa del nebulizador	≥ 0.2 ml/min
Líquido residual	<0.2 ml




PORTABLE COMPRESSOR NEBULIZER

Penetrate into Bronchus & Lung 2-5µm



MICRON NEBULIZATION FINE PARTICLES



3.6-4.0µm Median particle size MMAD	75% Percentage of effective particles
--	--

QUIET ENVIRONMENT SILENT ATOMIZATION

Low noise design, easy for children to sleep



Nebulizer is suitable for auxiliary treatment of most **respiratory diseases**



1. Nebulzier machine
2. Medicine cup
3. Air tube
4. Adult mask
5. Child mask
6. Mouthpiece
7. Elastic strap

