

Nebulizer kompresor Procyon menyediakan terapi inhalasi nebulisasi yang cepat, efisien, dan andal, dirancang untuk digunakan di rumah. Ini adalah nebulizer kompresor yang tenang dan bekerja keras, cocok untuk nebulisasi sebagian besar obat aerosol yang diresepkan untuk penyakit pernapasan umum. Nebulizer ini menawarkan aliran liter yang lebih tinggi, waktu perawatan yang lebih singkat, dan tingkat suara yang lebih rendah dibandingkan banyak nebulizer kompresor ekonomis lainnya.

**Fitur Kinerja dan Desain:**  
 Tekanan dan aliran yang andal  
 Nebulizer yang dapat digunakan kembali dapat disertakan  
 Jarak kaki kecil  
 Desain kecil dan kompak  
 Pegangan pembawa yang nyaman  
 Bagian mekanis logam

Weight: 2.18kg



Seri / Model No.	Seri Procyon / 70802
Dimensi Satuan	191mm*173mm*102mm
Berat	2.18kg
Persyaratan elektrik	100-240VAC, 50-60HZ
Tingkat Kebisingan Penyempresian	≤ 60dBA
Tekanan Maksimum	≥ 33psi
Aliran udara bebas kompresor	≥ 6L/menit
Ukuran partikel	2µm - 5µm
Tingkat Nebuliser	≥ 0.2ml/menit
Cairan Sisa	10.2ml




## PORTABLE COMPRESSOR NEBULIZER

Penetrate into Bronchus & Lung 2-5µm



## MICRON NEBULIZATION FINE PARTICLES



<b>3.6-4.0µm</b> Median particle size MMAD	<b>75%</b> Percentage of effective particles
--	--

## QUIET ENVIRONMENT SILENT ATOMIZATION

Low noise design, easy for children to sleep



Nebulizer is suitable for auxiliary treatment of most **respiratory diseases**



1. Nebulzier machine
2. Medicine cup
3. Air tube
4. Adult mask
5. Child mask
6. Mouthpiece
7. Elastic strap

