

Patient Turner 72110

SPECIFICATIONS

| | |
|-------------------------|---------------|
| Height | 40.9" - 49.6" |
| Leg Support Height Min | 6.3" |
| Leg Supportt Height Max | 7.1" |
| Base Diameter | 16.1" |
| Weight | 36lbs |
| Weight Capacity | 330lbs |
| Casters | 4" |
| Pad Dimensions | 6.9" * 4.9" |



330 lbs
weight capacity


DETAILED DESCRIPTION



The 72110 is a patient turner which supports the carer and service user to participate in short distance transfers, such as; seat to seat, bed to seat, wheelchair to toilet and so for. The patient turner allows the user to stand and turn whilst supporting them on short distance transfers.



The 72110 can adjust in height to help suit the needs of the user, whilst the handle has various gripping points for both the carer and user.





Two rear 8 cm (7.6 cm) wheels for effortless,
even with a heavy users - Easy to dismantle/stroller
attachment without tools for transportation and storage.



Padded lower leg supports can be adjusted in width, angle and height for user comfort. (tool free)



Low profile base facilitates user feet placement,
Substantial base area provides stability,
with rear bar also allowing carer to counterbalance.

SIMILAR PRODUCTS

