

# Power Hoyer Lift 71940



Weight Capacity  
400lbs / 181kg



Emergency Stop Button

Sturdy nylon sling

Rear casters with locking

## Product Description

Heavy duty construction of the Electric Patient Lift from Careage assists in a safe lift for both caregiver and patient.

This lift is ideal for home use. The wide base of the lift offers excellent stability while still being able to clear a low bed.

The battery is easy to charge and maintain by plugging the lift directly into the wall and not having to remove the battery from the device.

## Features:

1. Battery operated lift provides safe patient transfers
2. Emergency stop button and emergency lowering system
3. Heavy Duty steel construction patient lift supports up to 400 lbs.
4. Extra wide 6-point cradle with 360° rotation
5. Wide base provides additional stability, Manual lever opens and closes base legs
6. Base legs adjust easily and lock securely into open position with shifter handle

## Applications:



Lifting from floor



Bed-Chair transfer




Toileting/Bathing

## Detailed Description



Lift from floor, bed, commode, or wheelchair.

Easy to use Lift system has a six point cradle that will accommodate a six point sling, 4 point sling or a 2 point sling.





Easily widen and narrow the base by pulling the adjustment handle back to unlock.

The 400 lbs safe working load and greater spatial area make this lift a perfect fit for larger residents, fully meet the normal needs.



Remote control

ergonomic and easy to use

## Product Comparison:



	Battery Powered Patient Lift 71920	Battery Powered Patient Lift 71921	Bariatric Battery Powered Lift 71922	Electric Patient Lift 71940
 Weight Capacity	450 lbs/204kg	450 lbs/204kg	600 lbs/272kg	400 lbs/181kg
 Net Weight	49kg	47.6kg	59.9kg	39.9kg
 Gross weight	55.3kg	54.4kg	67.1kg	48.5kg
 Base Height	4.75"	4.5"	4.5"	4.4"
 Boom Height	27.5" - 76"	29" - 80.5"	32" - 75.5"	20" - 65"
 Casters	3"Front/5"Rear	3"Front/4.5"Rear	3" Front/4" Rear	4"
 Min and Max Base Width	41"/24"	41"/24"	45.75" / 25"	33" / 23"
 Weight Scale	Optional	/	Optional	/