

Electric Wheelchair Manufacturer

TRAVEL BUTLER SMART ELECTRIC WHEELCHAIR



PRODUCT STRUCTURE



Motor 1
Motor 2
500W dual motor,
DC motor /strong climbing ability



Carbonic acid battery
Driving distance is about 12km

PRODUCT DETAILS



Indicator light

Power switch

Horn


Power display

Speed adjustment

360 CONTROLLER

Wheelchair stop immediately when the hand leaves the controller.

PRODUCT ADVANTAGES



EABS Brake system

500W dual motor

High battery life

360°joystick



Hand/electric mode

Easy to fold

Carbon steel material

Anti-backover design

FOLDING DESIGN

Lightweight foldable portable design
Can be put in the trunk, can be boarded on the plane, can board the high-speed rail.



LEFT OR RIGHT HAND ADJUSTABLE

Can be adjusted as needed, the default is installed on the right



THE ARMREST CAN BE RAISED

Handrails can swing up and down
Easy to hid the patient up and down wheelchair



PRODUCT PARAMETER

Parameter	Category
Series/Model No.	Alphard Series/74502
Frame material	Carbon steel
Unfolded Dimension	1070MM*650MM*930MM
Folded Dimension	785mm*445mm*745mm
Motor	DC 24V、2*250W
Battery	DC 24V 12AH
Weight	36KG
Loading Capacity	120KG
Front wheel specifications	10" PU tire
Rear wheel specifications	16" Pneumatic tire
Max Speed	6KM/H
Driving Range	12KM

HAND ELECTRIC SWITCH

Easily switch between manual and electric, multiple ways to travel

