

# **Electric Wheelchair Manufacturer**







# TRAVEL BUTLER SMART ELECTRIC WHEELCHAIR



## PRODUCT STRUCTURE



Motor 1  
Motor 2  
500W dual motor,  
DC motor /strong climbing ability



Carbonic acid battery  
Driving distance is about 12km

## PRODUCT DETAILS



## PRODUCT ADVANTAGES



EABS Brake system




500W dual motor



High battery life



360°joystick



Hand/electric mode



Easy to fold



Carbon steel material



Anti-backover design

### FOLDING DESIGN

Lightweight foldable portable design  
Can be put in the trunk, can be boarded on the plane, can board the high-speed rail.



### LEFT OR RIGHT HAND ADJUSTABLE

Can be adjusted as needed, the default is installed on the right



### THE ARMREST CAN BE RAISED

Handrails can swing up and down  
Easy to hld the patient up and down wheelchair



### HAND ELECTRIC SWITCH

Easily switch between manual and electric, multiple ways to travel



## PRODUCT PARAMETER

Parameter	Category
Series/Model No.	Alphard Series/74502
Frame material	Carbon steel
Unfolded Dimension	1070MM*650MM*930MM
Folded Dimension	785mm*445mm*745mm
Motor	DC 24V、2*250W
Battery	DC 24V 12AH
Weight	36KG
Loading Capacity	120KG
Front wheel specifications	10" PU tire
Rear wheel specifications	16" Pneumatic tire
Max Speed	6KM/H
Driving Range	12KM