

Stair Climbing Wheelchair

Stair climbing wheelchair is mainly used for high-rise buildings to transfer patient up and down stairs. Suitable for home, hospitals, emergency centers etc. Foldable design and one person operate easily. Upgraded seat belt to ensure safety & anti-dropping. This stair lift is lightweight, non-broken, and quite labor-saving.



Stair Climbing Wheelchair Buying Guide



Weight Capacity
(265lbs)



Lithium battery



Tank track



Foldable



Product Structure

- | | |
|---------------------------|-----------------------------|
| (1) Smart operation panel | (7) Lithium battery |
| (2) Adjust the lever | (8) Climbing track |
| (3) Cushion/cushion | (9) Stretcher |
| (4) Rear handle | (10) Pedal |
| (5) Armrest | (11) Wheel |
| (6) Seat belt | (12) Anti-back wheel device |

Product Advantage

Transfer Easily

Safe and effectively transports patients up and downstairs and across the level ground as a wheelchair.

Easy to operate

Quickly and easily you can fold and store it in the car boot and for transportation.

Safe & Durable

Aluminum Alloy design make it durable and safe for transfer; The unique track structure is convenient and safe when the stretcher goes down the stairs. Dual-use as an electric wheelchair or emergency stretcher (Pull out the back rail).