

Máy phun khí dung nén Procyon cung cấp liệu pháp xông khí dung nhanh chóng, hiệu quả và đáng tin cậy, được thiết kế để sử dụng tại nhà. Đây là máy phun khí dung dung nén hoạt động êm ái và thích hợp để phun hầu hết các loại thuốc khí dung được kê đơn cho các bệnh về đường hô hấp thông thường. Nó cung cấp lưu lượng lớn cao hơn, thời gian sử lý ngắn hơn và mức âm thanh thấp hơn so với nhiều máy phun khí dung nén tiết kiệm khác.

Tính năng hiệu suất và thiết kế:

- Áp suất và dòng chảy đáng tin cậy
- Có thể kèm theo máy phun sương có thể tái sử dụng
- Điều chỉnh nhỏ
- Thiết kế nhỏ gọn
- Tay cầm tiện lợi
- Bộ phận cơ khí kim loại

Weight: 2.18kg



Sê-ri/Mẫu số	Dòng Procyon / 70802
Kích thước đơn vị	191mm*170mm*102mm
Cân nặng	2.18kg
Yêu cầu về điện	100-240VAC, 50-60HZ
Độ ồn khi vận hành	≤ 60 dBA
Áp suất tối đa	≤ 35psi
Lưu lượng khí từ đầu của máy nén	≤ 6L/phút
Kích thước hạt	2µm - 5µm
Tỷ lệ máy phun sương	≤ 0.2ml/phút
Chất lỏng dự	70.5ml




## PORTABLE COMPRESSOR NEBULIZER

Penetrate into Bronchus & Lung 2-5µm



Adults & Children

## MICRON NEBULIZATION FINE PARTICLES



<p><b>3.6-4.0µm</b> Median particle size MMAD</p>	<p><b>75%</b> Percentage of effective particles</p>
---	---

## QUIET ENVIRONMENT SILENT ATOMIZATION

Low noise design, easy for children to sleep



Nebulizer is suitable for auxiliary treatment of most **respiratory diseases**



1. Nebulzier machine
2. Medicine cup
3. Air tube
4. Adult mask
5. Child mask
6. Mouthpiece
7. Elastic strap

